

TERF2IP Polyclonal Antibody

Catalog No: #31399

Package Size: #31399-1 50ul #31399-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TERF2IP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TERF2IP (NP_061848.2).
Other Names	TERF2IP;DRIP5;RAP1
Accession No.	Uniprot:Q9NYB0GeneID:54386
Uniprot	Q9NYB0
GeneID	54386
Calculated MW	55kDa
SDS-PAGE MW	58KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of U2OS cells using TERF2IP antibody.



Western blot analysis of extracts of various cell lines, using TERF2IP antibody.

Background

The gene encodes a protein that is part of a complex involved in telomere length regulation. Pseudogenes are present on chromosomes 5 and 22.

Note: This product is for in vitro research use only