

RUFY2 Polyclonal Antibody

Catalog No: #31401



Package Size: #31401-1 50ul #31401-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RUFY2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RUFY2 (NP_001035882.1).
Other Names	RUFY2;RABIP4R;ZFYVE13
Accession No.	Uniprot:Q8WXA3GenelD:55680
Uniprot	Q8WXA3
GenelD	55680
Calculated MW	75kDa
SDS-PAGE MW	75kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RUFY2 antibody.

Note: This product is for in vitro research use only