

NLGN4X Polyclonal Antibody

Catalog No: #31403



Package Size: #31403-1 50ul #31403-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NLGN4X Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NLGN4X (NP_065793.1).
Other Names	NLGN4X; ASPGX2; AUTSX2; HLNX; HNL4X; NLGN4; neuroligin-4, X-linked
Accession No.	Swiss-Prot#:Q8N0W4NCBI Gene ID:57502
Uniprot	Q8N0W4
GeneID	57502;
Calculated MW	92kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NLGN4X at 1:1000 dilution.

Background

This gene encodes a member of the type-B carboxylesterase/lipase protein family. The encoded protein belongs to a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses. The encoded protein interacts with discs large homolog 4 (DLG4). Mutations in this gene have been associated with autism and Asperger syndrome. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only