

## ankyrin repeat domain 45 antibody

Catalog No: #23181

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	ankyrin repeat domain 45 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 58 and 266 of Human ANKRD45
Target Name	ankyrin repeat domain 45
Accession No.	Swiss-Prot:Q5TZF3Gene ID:339416
Uniprot	Q5TZF3
GeneID	339416;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 30kd

Western blotting: 1:500-1:3000

## Images



Sample(30 ug whole cell lysate)  
A: Hep G2  
12% SDS PAGE  
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only