

TNS3 Polyclonal Antibody

Catalog No: #31406



Package Size: #31406-1 50ul #31406-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TNS3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TNS3 (NP_073585.8).
Other Names	TNS3; TEM6; TENS1; tensin-3
Accession No.	Swiss-Prot#:Q68CZ2NCBI Gene ID:64759
Uniprot	Q68CZ2
GeneID	64759;
Calculated MW	190kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of A-549 cells, using TNS3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only