

NDUFS3 Polyclonal Antibody

Catalog No: #31419



Package Size: #31419-1 50ul #31419-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NDUFS3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NDUFS3 (NP_004542.1).
Other Names	NDUFS3;CI-30
Accession No.	Uniprot:O75489GeneID:4722
Uniprot	O75489
GeneID	4722
Calculated MW	30kDa
SDS-PAGE MW	30kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NDUFS3 antibody.

Background

This gene encodes one of the iron-sulfur protein (IP) components of mitochondrial NADH:ubiquinone oxidoreductase (complex I). Mutations in this gene are associated with Leigh syndrome resulting from mitochondrial complex I deficiency.

Note: This product is for in vitro research use only