

TRMT61B Polyclonal Antibody

Catalog No: #31429

Package Size: #31429-1 50ul #31429-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TRMT61B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TRMT61B (NP_060380.3).
Other Names	TRMT61B; tRNA methyltransferase 61B
Accession No.	Swiss-Prot#:Q9BVS5NCBI Gene ID:55006
Uniprot	Q9BVS5
GeneID	55006;
Calculated MW	52kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse liver, using TRMT61B at 1:1000 dilution.

Background

Note: This product is for in vitro research use only