

Gm13125 Polyclonal Antibody

Catalog No: #31436

Package Size: #31436-1 50ul #31436-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Gm13125 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse Gm13125 (NP_001108549.1).
Other Names	PRAMEF20; PRAMEF21; PRAME family member 20
Accession No.	Swiss-Prot#:Q5VT98NCBI Gene ID:627009
Uniprot	Q5VT98
GeneID	645425;
Calculated MW	55kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:100 - 1:200


Images



Immunohistochemistry of paraffin-embedded human liver cancer using Gm13125 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using Gm13125 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using Gm13125 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only