

FKBP2 Polyclonal Antibody

Catalog No: #31443



Package Size: #31443-1 50ul #31443-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FKBP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human FKBP2 (NP_004461.2).
Other Names	FKBP2;FKBP-13;PPIase
Accession No.	Uniprot:P26885GenelD:2286
Uniprot	P26885
GenelD	2286
Calculated MW	14kDa
SDS-PAGE MW	14kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000IHC 1:100 - 1:200

Images



Western blot analysis of extracts of various cell lines, using FKBP2 antibody.

Background

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It is thought to function as an ER chaperone and may also act as a component of membrane cytoskeletal scaffolds. Multiple alternatively

spliced variants, encoding the same protein, have been identified.

Note: This product is for in vitro research use only