

HOOK2 Polyclonal Antibody

Catalog No: #31470



Package Size: #31470-1 50ul #31470-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HOOK2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human HOOK2 (NP_001093646.1).
Other Names	HOOK2; HK2; protein Hook homolog 2
Accession No.	Swiss-Prot#:Q96ED9NCBI Gene ID:29911
Uniprot	Q96ED9
GeneID	29911;
Calculated MW	100kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using HOOK2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only