

NIN Polyclonal Antibody

Catalog No: #31472



Package Size: #31472-1 50ul #31472-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NIN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NIN (NP_891991.1).
Other Names	NIN;SCKL7;ninein
Accession No.	Uniprot:Q8N4C6GeneID:51199
Uniprot	Q8N4C6
GeneID	51199
Calculated MW	190kDa
SDS-PAGE MW	190kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000IHC 1:50 - 1:100


Images




Western blot analysis of extracts of Jurkat cells, using NIN antibody.



Immunohistochemistry of paraffin-embedded human stomach using NIN antibody.



Immunohistochemistry of paraffin-embedded mouse lung using NIN antibody.



Immunohistochemistry of paraffin-embedded rat kidney using NIN antibody.

Background

This gene encodes one of the proteins important for centrosomal function. This protein is important for positioning and anchoring the microtubules minus-ends in epithelial cells. Localization of this protein to the centrosome requires three leucine zippers in the central coiled-coil domain. Multiple alternatively spliced transcript variants that encode different isoforms have been reported.

Note: This product is for in vitro research use only