

# RNF31 Polyclonal Antibody

Catalog No: #31476



Package Size: #31476-1 50ul #31476-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	RNF31 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human RNF31 (NP_060469.4).
Other Names	RNF31; HOIP; Paul; ZIBRA; ring finger protein 31
Accession No.	Swiss-Prot#:Q96EP0NCBI Gene ID:55072
Uniprot	Q96EP0
GeneID	55072;
Calculated MW	125kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

## Images



Western blot analysis of extracts of various cell lines, using RNF31 at 1:1000 dilution.

## Background

The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. The encoded protein is the E3 ubiquitin-protein ligase component of the linear ubiquitin chain assembly complex. Two transcript variants encoding different isoforms have been found for this gene.

---

Note: This product is for in vitro research use only