

ITPK1 Polyclonal Antibody

Catalog No: #31504

Package Size: #31504-1 50ul #31504-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ITPK1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ITPK1 (NP_055031.2).
Other Names	ITPK1; ITRPK1; inositol-tetrakisphosphate 1-kinase
Accession No.	Swiss-Prot#:Q13572NCBI Gene ID:3705
Uniprot	Q13572
GeneID	3705;
Calculated MW	39kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ITPK1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only