

RELN Antibody

Catalog No: #46195

Package Size: #46195-1 50ul #46195-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RELN Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Rt
Specificity	RELN Antibody detects endogenous levels of total RELN
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RELN
Target Name	RELN
Other Names	RELN, LIS2, PRO1598, Reelin, RL
Accession No.	Swiss-Prot#: P78509NCBI Gene ID:5649
Uniprot	P78509
GeneID	5649;
Calculated MW	388kDa
Concentration	1.0mg/mL
Formulation	PBS with 0.05% proclin300, 50% glycerol, pH7.3.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of extracts of HepG2 cells, using RELN Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Exposure time: 3min.

Note: This product is for in vitro research use only