

SFI1 Polyclonal Antibody

Catalog No: #31522



Package Size: #31522-1 50ul #31522-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SFI1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SFI1 (NP_055590.2).
Other Names	SFI1; PISD; PPP1R139; hSfi1p; protein SFI1 homolog
Accession No.	Swiss-Prot#:A8K8P3NCBI Gene ID:9814
Uniprot	A8K8P3
GeneID	9814;
Calculated MW	148kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of MCF-7 cells, using SFI1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only