

PRDM11 Polyclonal Antibody

Catalog No: #31576



Package Size: #31576-1 50ul #31576-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PRDM11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PRDM11 (NP_001243625.1).
Other Names	PRDM11; PFM8; PR/SET domain 11
Accession No.	Swiss-Prot#:Q9NQV5NCBI Gene ID:56981
Uniprot	Q9NQV5
GeneID	56981;
Calculated MW	58kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PRDM11 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only