

PPP4R4 Polyclonal Antibody

Catalog No: #31578

Package Size: #31578-1 50ul #31578-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PPP4R4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PPP4R4 (NP_478144.1).
Other Names	PPP4R4; CFAP14; KIAA1622; PP4R4; protein phosphatase 4 regulatory subunit 4
Accession No.	Swiss-Prot#:Q6NUP7NCBI Gene ID:57718
Uniprot	Q6NUP7
GeneID	57718;
Calculated MW	99kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse lung, using PPP4R4 at 1:1000 dilution.

Background

The protein encoded by this gene is a HEAT-like repeat-containing protein. The HEAT repeat is a tandemly repeated, 37-47 amino acid long module occurring in a number of cytoplasmic proteins. Arrays of HEAT repeats form a rod-like helical structure and appear to function as protein-protein interaction surfaces. The repeat-containing region of this protein has some similarity to the constant regulatory domain of the protein phosphatase 2A PR65/A subunit. The encoded protein binds protein serine/threonine phosphatase 4c in the cytoplasm.

Note: This product is for in vitro research use only