

C19orf40 Polyclonal Antibody

Catalog No: #31586



Package Size: #31586-1 50ul #31586-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	C19orf40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human C19orf40 (NP_689479.1).
Other Names	FAAP24; C19orf40; Fanconi anemia core complex-associated protein 24
Accession No.	Swiss-Prot#:Q9BTP7NCBI Gene ID:91442
Uniprot	Q9BTP7
GeneID	91442;
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of MCF-7 cells, using C19orf40 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only