

METTL14 Polyclonal Antibody

Catalog No: #31591



Package Size: #31591-1 50ul #31591-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	METTL14 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human METTL14 (NP_066012.1).
Other Names	METTL14; hMETTL14; methyltransferase like 14
Accession No.	Swiss-Prot#:Q9HCE5NCBI Gene ID:57721
Uniprot	Q9HCE5
GeneID	57721;
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using METTL14 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only