

RHG28 Antibody

Catalog No: #46205



Package Size: #46205-1 50ul #46205-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RHG28 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	RHG28 Antibody detects endogenous levels of total RHG28
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RHG28
Target Name	RHG28
Other Names	ARHGAP28, KIAA1314
Accession No.	Swiss-Prot#: Q9P2N2
Uniprot	Q9P2N2
Calculated MW	82 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only