

HDGFRP3 Polyclonal Antibody

Catalog No: #31640



Package Size: #31640-1 50ul #31640-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HDGFRP3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human HDGFRP3 (NP_057157.1).
Other Names	HDGFRP3; CGI-142; HDGF-2; HDGF2; HRP-3; hepatoma-derived growth factor-related protein 3
Accession No.	Swiss-Prot#:_NCBI Gene ID:50810
GenelD	:50810
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using HDGFRP3 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat spinal cord using HDGFRP3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human uterine cancer using HDGFRP3 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only