

EXOSC4 Polyclonal Antibody

Catalog No: #31675



Package Size: #31675-1 50ul #31675-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EXOSC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human EXOSC4 (NP_061910.1).
Other Names	EXOSC4; RRP41; RRP41A; Rrp41p; SKI6; Ski6p; hRrp41p; p12A; exosome component 4
Accession No.	Swiss-Prot#:Q9NPD3NCBI Gene ID:54512
Uniprot	Q9NPD3
GeneID	54512;
Calculated MW	26kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:100 - 1:200

Images



Western blot analysis of extracts of various cell lines, using EXOSC4 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only