

FRMD6 Polyclonal Antibody

Catalog No: #31802



Package Size: #31802-1 50ul #31802-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FRMD6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human FRMD6 (NP_689543.1).
Other Names	FRMD6; C14orf31; EX1; Willin; c14_5320; FERM domain containing 6
Accession No.	Swiss-Prot#:Q96NE9NCBI Gene ID:122786
Uniprot	Q96NE9
GeneID	122786;
Calculated MW	72kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:1000 - 1:4000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using FRMD6 at 1:1000 dilution.




Immunofluorescence analysis of C6 cells using FRMD6 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of U2OS cells using FRMD6 at dilution of 1:100. Blue: DAPI for nuclear staining.




Immunohistochemistry of paraffin-embedded Rat heart using FRMD6 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human thyroid cancer using FRMD6 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse heart using FRMD6 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using FRMD6 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only