

## SH3GL3 Antibody

Catalog No: #46239

Package Size: #46239-1 50ul #46239-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                        |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | SH3GL3 Antibody                                                                                                                        |
| Host Species          | Rabbit                                                                                                                                 |
| Clonality             | Polyclonal                                                                                                                             |
| Purification          | Immunogen affinity purified                                                                                                            |
| Applications          | WB                                                                                                                                     |
| Species Reactivity    | Hu Ms Rt                                                                                                                               |
| Specificity           | SH3GL3 Antibody detects endogenous levels of total SH3GL3                                                                              |
| Immunogen Type        | Peptide                                                                                                                                |
| Immunogen Description | A synthesized peptide derived from human SH3GL3                                                                                        |
| Target Name           | SH3GL3                                                                                                                                 |
| Other Names           | SH3GL3, CNSA3, EEN-2B-L3, EEN-B2, Endophilin-3, Endophilin-A3, HsT19371, SH3 domain protein 2C, SH3P13, SH3-domain GRB2-like 3, SH3D2C |
| Accession No.         | Swiss-Prot#: Q99963NCBI Gene ID:6457                                                                                                   |
| Uniprot               | Q99963                                                                                                                                 |
| GeneID                | 6457;                                                                                                                                  |
| Calculated MW         | 39 kDa                                                                                                                                 |
| Concentration         | 1.0mg/mL                                                                                                                               |
| Formulation           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                         |
| Storage               | Store at -20°C                                                                                                                         |

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 40% over night

---

Note: This product is for in vitro research use only