

SLC9A2 Antibody

Catalog No: #46256

Package Size: #46256-1 50ul #46256-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC9A2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SLC9A2 Antibody detects endogenous levels of total SLC9A2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SLC9A2
Target Name	SLC9A2
Other Names	SLC9A2, Na(+)/H(+) exchanger 2, NHE-2, Sodium-hydrogen exchanger 2, Sodium/hydrogen exchanger 2, NHE2
Accession No.	Swiss-Prot#: Q9UBY0NCBI Gene ID:6549
Uniprot	Q9UBY0
GeneID	6549;
Calculated MW	92 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only