

## TCEAL7 Antibody

Catalog No: #46283

Package Size: #46283-1 50ul #46283-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TCEAL7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	TCEAL7 Antibody detects endogenous levels of total TCEAL7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TCEAL7
Target Name	TCEAL7
Other Names	Transcription elongation factor A protein-like 7 TCEA-like protein 7 Transcription elongation factor S-II protein-like 7
Accession No.	Swiss-Prot#: Q9BRU2NCBI Gene ID:56849
Uniprot	Q9BRU2
GeneID	56849;
Calculated MW	12 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

WB 1:500-2000

## Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4o½ over night

---

Note: This product is for in vitro research use only