

SAMHD1 (Phospho-Thr592) Antibody

Catalog No: #12805



Package Size: #12805-1 50ul #12805-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SAMHD1 (Phospho-Thr592) Antibody
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	SAMHD1 (Phospho-Thr592) Antibody detects endogenous levels of SAMHD1 only when phosphorylated at Thr592
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human SAMHD1 around the phosphorylation site of Thr592.
Target Name	SAMHD1 (Phospho-Thr592)
Modification	Phospho
Other Names	SAMHD1, AGS5, CHBL2, DCIP, HDDC1, Monocyte protein 5, MOP-5, MOP5, SBBI88, Mg11, SAM domain and HD domain 1
Accession No.	Swiss-Prot#: Q9Y3Z3NCBI Gene ID: 25939
Uniprot	Q9Y3Z3
GeneID	25939;
Target Species	human
Calculated MW	72KD
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only