

LRP6 (Phospho-Ser1490) Antibody

Catalog No: #12809

Package Size: #12809-1 50ul #12809-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP6 (Phospho-Ser1490) Antibody
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	LRP6 (Phospho-Ser1490) Antibody detects endogenous levels of LRP6 only when phosphorylated at Ser1490
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human LRP6 around the phosphorylation site of Ser1490.
Target Name	LRP6 (Phospho-Ser1490)
Modification	Phospho
Other Names	LRP6, ADCAD2, LRP-6
Accession No.	Swiss-Prot#: O75581NCBI Gene ID: 4040
Uniprot	O75581
GeneID	4040;
Target Species	human
Calculated MW	180 KD
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

Note: This product is for in vitro research use only