

## BLBP Polyclonal Antibody

Catalog No: #46784

Package Size: #46784-1 50ul #46784-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BLBP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Hu Rt
Specificity	The antibody detects endogenous BLBP
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 80-132
Target Name	BLBP
Other Names	Fatty acid-binding protein, brain (Brain lipid-binding protein) (BLBP) (Brain-type fatty acid-binding protein) (B-FABP) (Fatty acid-binding protein 7) (Mammary-derived growth inhibitor related)
Accession No.	Swiss-Prot:O15540NCBI Gene ID:2173
Uniprot	O15540
GeneID	2173;
Calculated MW	15 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1B£B'1500-1B£B'12000

## Images



Western blot analysis of HEPG2 293T SHY5Y Jurkat lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000

---

Note: This product is for in vitro research use only