

## CD108 Polyclonal Antibody

Catalog No: #46809

Package Size: #46809-1 50ul #46809-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CD108 Polyclonal Antibody                                                                                                                                       |
| Host Species          | Rabbit                                                                                                                                                          |
| Clonality             | Polyclonal                                                                                                                                                      |
| Purification          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                       |
| Applications          | IHC                                                                                                                                                             |
| Species Reactivity    | Hu Ms Rt                                                                                                                                                        |
| Specificity           | The antibody detects endogenous CD108                                                                                                                           |
| Immunogen Type        | protein                                                                                                                                                         |
| Immunogen Description | Synthetic peptide from human protein at AA range: 101-150                                                                                                       |
| Target Name           | CD108                                                                                                                                                           |
| Other Names           | Semaphorin-7A (CDw108) (JMH blood group antigen) (John-Milton-Hargen human blood group Ag) (Semaphorin-K1) (Sema K1) (Semaphorin-L) (Sema L) (CD antigen CD108) |
| Accession No.         | Swiss-Prot:O75326NCBI Gene ID:8482                                                                                                                              |
| Uniprot               | O75326                                                                                                                                                          |
| GeneID                | 8482;                                                                                                                                                           |
| Concentration         | 1 mg/ml                                                                                                                                                         |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                         |
| Storage               | Store at -20°C                                                                                                                                                  |

## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200

---

Note: This product is for in vitro research use only