

Jun B Polyclonal Antibody

Catalog No: #46903



Package Size: #46903-1 50ul #46903-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Jun B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Jun B Polyclonal Antibody detects endogenous levels of Jun B
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human Jun B. at AA range: 231-280
Target Name	Jun B
Other Names	Transcription factor jun-B
Accession No.	Swiss-Prot:P17275NCBI Gene ID:3726
Uniprot	P17275
GeneID	3726;
Calculated MW	40 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of 1,mouse-kidney 2,mouse-heart cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B1£C, 1 hour)

Note: This product is for in vitro research use only