

AVPR1A Polyclonal Antibody

Catalog No: #46914



Package Size: #46914-1 50ul #46914-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	AVPR1A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	AVPR1A Polyclonal Antibody detects endogenous levels of AVPR1A
Immunogen Type	protein
Immunogen Description	Synthesized peptide derived from human AVPR1A. at AA range: 11-60
Target Name	AVPR1A
Other Names	Vasopressin V1a receptor (V1aR) (AVPR V1a) (Antidiuretic hormone receptor 1a) (Vascular/hepatic-type arginine vasopressin receptor)
Accession No.	Swiss-Prot:P37288NCBI Gene ID:552
Uniprot	P37288
GeneID	552;
Calculated MW	55 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B1500-1B£B12000

Images



Western Blot analysis of mouse-lung cells using primary antibody diluted at 1:500(°C overnight). Secondary antibodyB£B1Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25B'£C, 1 hour)

Note: This product is for in vitro research use only