

BTF3L4 Antibody

Catalog No: #46358

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	BTF3L4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total BTF3L4 protein.
Immunogen Type	peptide
Immunogen Description	Full length fusion protein of human BTF3L4
Target Name	BTF3L4
Accession No.	Swiss-Prot:Q96K17NCBI Gene ID:91408NCBI Protein:BC022371
Uniprot	Q96K17
GeneID	91408;
Calculated MW	17 kDa
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000


Immunohistochemistry: 1: 25-100

Images




Gel: 12%SDS-PAGE

lysate: 40 B₂g, Lane: Hela cell lysate,Primary antibody: 46358B£B`BTF3L4 Antibody) at dilution
1/350 dilution,Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 1 minute



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46358(BTF3L4 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 46358(BTF3L4 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only