

PGK1 (Phospho-Tyr324) Antibody

Catalog No: #13437

Package Size: #13437-1 50ul #13437-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PGK1 (Phospho-Tyr324) Antibody
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Applications	WB IHC
Species Reactivity	Human, Mouse, Rat
Specificity	Phospho-PGK1(Y324) Antibody detects endogenous levels of PGK1 only when phosphorylated at Y324.
Immunogen Description	A synthesized peptide derived from human PGK1 around the phosphorylation site of Y324.
Target Name	PGK1
Other Names	Cell migration-inducing gene 10 protein; Epididymis secretory sperm binding protein Li 68p; HEL S 68p; MGC117307; MGC8947; MIG10; pgk1; PGK1_HUMAN; PGKA; Phosphoglycerate kinase 1; Primer recognition protein 2; PRP 2;
Accession No.	Swiss-Prot#:P00558NCBI Gene ID:5230
Uniprot	P00558
GeneID	5230;
Calculated MW	45kda
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt.

Application Details

WB: 1:500-1:2000

IHC: 1:50-1:200

Note: This product is for in vitro research use only