

Human Coagulation Factor VIII (F8) ELISA Kit

Catalog No: #EK20106

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Coagulation Factor VIII (F8) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Hu
Other Names	Coagulation factor VIII,Antihemophilic factor,AHF,Procoagulant component,F8,F8C
Accession No.	UniProtKB: P00451 Gene ID:2157
Uniprot	P00451
GeneID	2157;
Storage	12 months

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.055ng/mL
Sample Type:tissue homogenates, cell lysates or other biological fluids.
Sample Volume: 100 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only