

ELP3 Polyclonal Antibody

Catalog No: #33183



Package Size: #33183-1 50ul #33183-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ELP3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ELP3 (NP_060561.3).
Other Names	ELP3; KAT9; elongator complex protein 3
Accession No.	Swiss-Prot#:Q9H9T3NCBI Gene ID:55140
Uniprot	Q9H9T3
GeneID	55140;
Calculated MW	62kDa
Storage	Store at -20oQ½? Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.3.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images



Western blot analysis of extracts of Mouse kidney, using ELP3 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat cerebellum using ELP3 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse cerebellum using ELP3 at dilution of 1:100 (40x lens).

Note: This product is for in vitro research use only