

CDRT15L2 Antibody

Catalog No: #46459

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CDRT15L2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CDRT15L2 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human CDRT15L2
Target Name	CDRT15L2
Accession No.	Swiss-Prot:A8MXV6 NCBI Gene ID:256223NCBI Protein:NP_001177719
Uniprot	A8MXV6
GeneID	256223;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 25-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 46459(CDRT15L2 Antibody) at dilution 1/25, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only