

Goat Anti-Rabbit IgG (H&L) Biotin

Catalog No: #L35041



Package Size: #L35041-1 100 ul #L35041-2 500 ul #L35041-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Goat Anti-Rabbit IgG (H&L) Biotin |
| Brief Description | Conjugated Secondary Antibody |
| Host Species | Goat |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | E, WB, IH, IC |
| Species Reactivity | Rb |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Rabbit IgG. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Rabbit IgG |
| Conjugates | Biotin |
| Formulation | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide |

Application Details

E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)

Note: This product is for in vitro research use only