

Rabbit Anti-Human IgG+IgM+IgA FITC

Catalog No: #L35086

Package Size: #L35086-1 100 ul #L35086-2 500 ul #L35086-3 1 ml

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Rabbit Anti-Human IgG+IgM+IgA FITC |
| Brief Description | Conjugated Secondary Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | IF, FC |
| Species Reactivity | H |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG IgM IgA. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG+IgM+IgA |
| Conjugates | FITC |
| Formulation | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide |

Application Details

IF (1/250 - 1/2000), FC (1/250 - 1/2000)

Note: This product is for in vitro research use only