

Mouse Anti-Human IgG Fab - AF790

Catalog No: #L35646



Package Size: #L35646-1 100 ul #L35646-2 500 ul #L35646-3 1 ml

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Mouse Anti-Human IgG Fab - AF790 |
| Brief Description | Conjugated Secondary Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Purification | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Applications | E, IH, IF, FC |
| Species Reactivity | H |
| Specificity | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fab. No antibody was detected against non immunoglobulin serum proteins. |
| Immunogen Description | Human IgG Fab |
| Conjugates | AF790 |
| Formulation | 1 mg/ml, liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |

Application Details

IH (1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), E (Use at an assay dependent concentration)

Note: This product is for in vitro research use only