

## HEATR9 Antibody

Catalog No: #46572

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HEATR9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total HEATR9 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the N terminal of human HEATR9
Target Name	HEATR9
Other Names	C17orf66
Accession No.	Swiss-Prot:A2RTY3NCBI Gene ID:256957NCBI Protein:BC033734
Uniprot	A2RTY3
GeneID	256957;
Concentration	2.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1: 50-300

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 46572(HEATR9 Antibody) at dilution 1/110, on the right is treated with fusion protein. (Original magnification: x200)

## Background

The function of this protein remains unknown.

---

Note: This product is for in vitro research use only