

## ZNF468 Antibody

Catalog No: #46721

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                 |
|-----------------------|-----------------------------------------------------------------|
| Product Name          | ZNF468 Antibody                                                 |
| Host Species          | Rabbit                                                          |
| Clonality             | Polyclonal                                                      |
| Purification          | Antigen affinity purification                                   |
| Applications          | IHC                                                             |
| Species Reactivity    | Hu                                                              |
| Specificity           | The antibody detects endogenous levels of total ZNF468 protein. |
| Immunogen Type        | peptide                                                         |
| Immunogen Description | Full length fusion protein of human ZNF468                      |
| Target Name           | ZNF468                                                          |
| Accession No.         | Swiss-Prot:Q5VIY5NCBI Gene ID:90333NCBI Protein:BC063143        |
| Uniprot               | Q5VIY5                                                          |
| GeneID                | 90333;                                                          |
| Concentration         | 0.8mg/ml                                                        |
| Formulation           | Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol. |
| Storage               | Store at -20°C                                                  |

## Application Details

Immunohistochemistry: 1: 30-150

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 46721(ZNF468 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x200)

## Background

May be involved in transcriptional regulation.

---

Note: This product is for in vitro research use only