

ZSCAN12 Antibody

Catalog No: #46723

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZSCAN12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZSCAN12 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic protein corresponding to residues near the C terminal of human ZSCAN12
Target Name	ZSCAN12
Other Names	ZFP96; ZNF96; ZNF305; ZNF29K1; dJ29K1.2
Accession No.	Swiss-Prot:O43309NCBI Gene ID:9753NCBI Protein:BC041661
Uniprot	O43309
GeneID	9753;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 40-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 46723(ZSCAN12 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using 46723(ZSCAN12 Antibody) at dilution 1/50, on the right is treated with fusion protein. (Original magnification: x200)

Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only