

FN1 Polyclonal Conjugated Antibody

Catalog No: #C40932

Package Size: #C40932-AF350 100ul #C40932-AF405 100ul #C40932-AF488 100ul

#C40932-AF555 100ul #C40932-AF594 100ul #C40932-AF647 100ul

#C40932-AF680 100ul #C40932-AF750 100ul #C40932-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FN1 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	FN1 Polyclonal Antibody detects endogenous levels of FN1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FN1.
Target Name	FN1
Other Names	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
Accession No.	Swiss-Prot: P02751NCBI Gene ID: 2335
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
SDS-PAGE MW	260kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Note: This product is for in vitro research use only