

MFAP3L (Phospho-Tyr287) Antibody

Catalog No: #SAB256P

Package Size: #SAB256P-1 50ul #SAB256P-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MFAP3L (Phospho-Tyr287) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	Custom antibody
Species Reactivity	Hu
Immunogen Type	Peptide-KLH
Target Name	MFAP3L
Other Names	Microfibrillar-associated protein 3-like; Testis development protein NYD-SP9
Accession No.	uniprot:O75121
Calculated MW	45kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Background

May participate in the nuclear signaling of EGFR and MAPK1/ERK2. May have a role in metastasis.

Note: This product is for in vitro research use only