

PDZ-GEF Antibody

Catalog No: #SAB404



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PDZ-GEF Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Fruit fly
Immunogen Type	Peptide-KLH
Target Name	PDZ-GEF
Other Names	D-PDZGEFoΩ'½oΩ'½DGEFoΩ'½oΩ'½dizzyoΩ'½oΩ'½Dmel\CG9491oΩ'½oΩ'½PDF-GEF Dizzy
Accession No.	uniprot:Q9VMF3
Calculated MW	172kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 69% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only