

ORF2 Antibody

Catalog No: #SAB442

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ORF2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Cherry mottle leaf virus
Immunogen Type	Peptide-KLH
Target Name	ORF2
Other Names	47kDa protein, ORF2
Accession No.	uniprot:Q9IMS8
Calculated MW	47kDa
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 90% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Western blotting: 1:500~1:1000

Note: This product is for in vitro research use only