

CX3CL1 Rabbit Polyclonal Antibody

Catalog No: #29282



Package Size: #29282-1 50ul #29282-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CX3CL1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB IHC IF
Species Reactivity	Human,Mouse,Rat
Other Names	CX3CL1;ABCD-3;C3Xkine;CXC3;CXC3C;NTN;NTT;SCYD1;fractalkine
Accession No.	Swiss-Prot#:P78423NCBI Gene ID:6376
Uniprot	P78423
GeneID	6376
Calculated MW	42kDa
Concentration	1.08mg/mL
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20oΩ%oΩ½. Avoid freeze / thaw cycles.


Application Details

WB	1:200 - 1:2000
IHC	1:50 - 1:200
IF	1:50 - 1:200


Images



Western blot analysis of extracts of rat brain, using CX3CL1 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded human brain using CX3CL1 Rabbit pAb (29282) at dilution of 1:150.



Immunofluorescence analysis of NIH-3T3 cells using CX3CL1 Rabbit pAb (29282) at dilution of 1:100 . Blue: DAPI for nuclear staining.

Note: This product is for in vitro research use only