

Recombinant human Probable histidine--tRNA ligase, mitochondrial

Catalog No: #AP70378

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP70378-1 20ug #AP70378-2 100ug #AP70378-3 1mg

Description

Product Name	Recombinant human Probable histidine--tRNA ligase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:34-506aaSequence Info:Full Length
Other Names	Histidine--tRNA ligase-likeHistidyl-tRNA synthetase ;HisRS
Accession No.	P49590
Uniprot	P49590
GeneID	23438;
Calculated MW	69.4 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	QSQVAEAVLTSQLKAHQEKPNIKTPKGTRDLSPQHMVVREKILDLVISCFKRHGAKGMDTPAFELKETLTEK YGEDSGLMYDLKDQGGELLSLRYDLTVPFARYLAMNKVKMKRYHVGKVWRRESPTIVQGRYREFCQCDFD IAGQFDPMIPDAECLKIMCEILSGLQLGDFLIKVNDRRIVDGMFAVCGVPESKFRAICSSIDKLDKMAWKDVRHE MVKKGLAPEVADRIGDYVQCHGGVSLVEQMFQDPRLSQNKQALEGLGDLKLLFEYLTLFGIADKISFDLSLA RGLDYTGVIYEAVLLQTPTQAGEEPLNVGSVAAGGRYDGLVGMFDPKGHKVPCVGLSIGVERIFYIVEQRMK TKGEKVRTTETQVFVATPQKNFLQERLKLIAELWDSGIKAEMLYKNNPKLLTQLHYCESTGIPLVVIIGEQLKE GVIKIRSVASREEVAIKRENFVAEIQRLSES
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

An enzyme assisted RP-RPLC approach for in-depth analysis of human liver phosphoproteome.Bian Y., Song C., Cheng K., Dong M., Wang F., Huang J., Sun D., Wang L., Ye M., Zou H.J. Proteomics 96:253-262(2014)Research Topic:Metabolism

Note: This product is for in vitro research use only