

Recombinant mouse Importin subunit beta-1

Catalog No: #AP70480



Package Size: #AP70480-1 20ug #AP70480-2 100ug #AP70480-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant mouse Importin subunit beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-876aaSequence Info:Full Length
Other Names	Karyopherin subunit beta-1Nuclear factor p97Pore targeting complex 97KDA subunit ;PTAC97SCG
Accession No.	P70168
Uniprot	P70168
GeneID	16211;
Calculated MW	101.2 kDa
Tag Info	N-terminal 6xHis-tagged
Target Sequence	MELITILEKTVSPDRLELEAAQKFLERAAVENLPTFLVELSRVLANPGNSQVARVAAGLQIKNSLTSKDPDIKAQ YQQRWLAI DANARREVKNYVLQTLGTETYPSSASQCVAGIACAEIPVSQWPELIPQLVANVTNPNSTEHMKE STLEAIGYICQDIDPEQLQDKSNEILTAIIQGMRKEEPSNNVKLAATNALLNSLEFTKANFDKESERHFIMQVCE ATQCPDTRVRVAALQNLVKIMSLYYQYMETYMGPALFAITIEAMKSDIDEVALQGIEFWSNVCDEEMDLAIEAS EAAEQGRPPEHTSKFYAKGALQYLVPILTQTLTKQDENDDDDWNPCKAAGVCLMLLSTCCEDDIVPHVLPFI KEHIKNPDWRYRDAAVMAFGSILEGPEPNQLKPLVIQAMPTLIELMKDPSVVVRDTTAWTVGRICELLPEAAIN DVYLAPLLQCLIEGLSAEPRVASNVCWAFSSLAEEAAYEADVADDQEEPATYCLSSSFELIVQKLETTDRPDG HQNNLRSSAYESLMEIVKNSAKDCYPVQKTTLVIMERLQQVLQMESHQSTSDRIQFNDLQSLLCATLQNVL KVQHQDALQISDVVMASLLRMFQSTAGSGGVQEDALMAVSTLVEVLGGEFLKYMEAFKPFGLGILKNYAEYQ VCLAAVGLVGDLCRALQSNILPFCDEVMQLLENLGNENVHRSVKPQILSVFGDIALAIGGEFKKYLEVVLNTLQ QASQAQVDKSDFDMDYLNELRESCLEAYTGIVQGLKGDQENVHPDVMLVQPRVEFILSFIDHIAGDEDHTDG VVACAAGLIGDLCTAFGKDVLKLVEARPMIHELLTEGRRSKTNKAKTLATWATKELRKLKNQA
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Functions in nuclear protein import, either in association with an adapter protein, like an importin-alpha subunit, which binds to nuclear localization signals (NLS) in cargo substrates, or by acting as autonomous nuclear transport receptor. Acting autonomously, serves itself as NLS receptor. Docking of the importin,substrate complex to the nuclear pore complex (NPC) is mediated by KPNB1 through binding to nucleoporin FxFG repeats and the complex is subsequently translocated through the pore by an energy requiring, Ran-dependent mechanism. At the nucleoplasmic side of the NPC, Ran binds to importin-beta and the three components separate and importin-alpha and -beta are re-exported from the nucleus to the cytoplasm where GTP hydrolysis releases Ran from importin. The directionality of nuclear import is thought to be conferred by an asymmetric distribution of the

GTP- and GDP-bound forms of Ran between the cytoplasm and nucleus. Mediates autonomously the nuclear import of ribosomal proteins RPL23A, RPS7 and RPL5. Binds to a beta-like import receptor binding (BIB) domain of RPL23A. In association with IPO7 mediates the nuclear import of H1 histone. In vitro, mediates nuclear import of H2A, H2B, H3 and H4 histones. In case of HIV-1 infection, binds and mediates the nuclear import of HIV-1 Rev. Imports SNAI1 and PRKCI into the nucleus .

References

Novel importin-alpha family member Kpna7 is required for normal fertility and fecundity in the mouse. Hu J., Wang F., Yuan Y., Zhu X., Wang Y., Zhang Y., Kou Z., Wang S., Gao S.J. Biol. Chem. 285:33113-33122(2010) Research Topic: Others

Note: This product is for in vitro research use only